

## Comparison of several approaches for vitamin D metabolite analysis

Xuejun (June) Zang, Chhavi Bhardwaj, Slobodan Milasinovic, Anil Oroskar, Asha Oroskar  
**Orochem Technologies, 340 Shuman Blvd., Naperville, IL 60563**

### Introduction


Vitamin D is a fat soluble vitamin responsible for regulating many physiological functions, among them, the absorption of calcium and phosphate. While the importance of Vitamin D has been recognized in maintaining the overall health and well being of individuals, increased clinical awareness of the prevalence of vitamin D deficiency and insufficiency has led to a need for large volume vitamin D metabolite assay.

Vitamin D and its metabolites circulate in blood bound to serum proteins. Thus to be able to analyze the metabolites, one will need to dissociate them from the serum proteins. LC-MS/MS has emerged as a reliable method of choice for vitamin D metabolite analysis due to its high specificity and high accuracy at low ng/mL level. However, LC-MS/MS is easily affected by phospholipids content in the serum samples.

We present two different approaches here for efficient vitamin D analysis:

- 1 Vitamin D metabolites extraction plate, It is simple and fast, while provides reasonable recoveries.
- 2 Sagacity HL SPE products, provides high recoveries and reproducibility of mono- and di- hydroxyl vitamin D.

### Oroflex Robotic Workflow



### Experiments

**Instruments** All automated extractions were carried out using Orochem's Oroflex Personal Pipettor robot. All LC-MS/MS method used a Shimadzu HPLC system coupled to an API 3000 mass spectrometer with a turbo spray ESI source operated in positive mode.

**Materials** Orochem's 96-well Vitamin D Metabolite Extraction plate, and Sagacity HL SPE plate were used for all extractions. The 25-hydroxyvitamin D2 and D3, 1,25-dihydroxyvitamin D3 standards were purchased from Sigma-aldrich. An EZYPRESS HT 96-well plate positive pressure manifold unit was used for conditioning, washing and processing the Vitamin D metabolite extraction plate and Sagacity HL SPE plate. Vitamin D and metabolites free human serum was purchased from Golden West Biologicals (Temecula, CA)

#### Procedures

##### Vitamin D metabolites extraction plate

1. Aliquot extraction reagent into each well
2. Add 0.1 mL of fortified serum sample
3. Wait about 1-2 minutes, then apply positive pressure to the plate such that fluid flow through the extraction plate drop by drop, total process is less than 10 minutes for 96-well plate
4. The vitamin D eluate is then evaporated and reconstituted with 0.1 mL of a solution of 80% methanol in water

##### Sagacity HL SPE plate

1. Condition the plate with 1 mL of acetonitrile followed by 1m of water.
2. Pretreat the fortified serum 0.1 mL with 50  $\mu$ L of isopropyl alcohol first, then add 70  $\mu$ L of water.
3. Load pretreated serum to the SPE plate, and wash with 1 mL of water followed by 1mL of 30% acetonitrile in water.
4. The metabolites are then eluted with 1 mL of acetonitrile. The vitamin D eluate was evaporated and reconstituted with 0.1 mL solution of 50% methanol in water.

**LC-MS/MS analysis** used Orochem's Reliasil C18 column, 3  $\mu$ m, 50 x 4.6 mm with a 87% Methanol (0.1% Formic acid, 15 mM Ammonium Acetate) mobile phase.

### Results

#### 1. Extraction Reagent and Recovery of Vitamin D Metabolites

We evaluate variations of extraction reagents (methanol, acetonitrile, different ratio of methanol and acetonitrile, with or without formic acid) and different ratio of reagent to serum (3:1, 4:1). The results show that methanol/acetonitrile 1:1 and reagent to serum 4:1 ratio provides the best recovery and phospholipids removal rate.

#### 2. Comparison of Vitamin D Metabolite Extraction Plate with Commercial Phospholipid Removal Plate


We evaluate three commercial phospholipid removal plates (B-D) with Vitamin D Metabolites Extraction Plate (A). For our plate, we use 4:1 reagent mixture (methanol/acetonitrile, 1:1) to serum ratio. We followed the instructions of each other commercial plates for its extraction.

**Table 1** - Comparison of Vitamin D Recovery for Commercial Kits

Kit	1,25-dihydroxy vitamin D3 Recovery	25-hydroxy vitamin D2 Recovery
A	74.0%	73.0%
B	36.9%	14.0%
C	52.8%	40.5%
D	63.5%	67%

Compared to some commercial phospholipid depletion plate, vitamin D metabolites gave the best recoveries.

**Figure 1** - Phospholipid profile of serum extracts.



#### 3. Sagacity HL SPE Plate

When using Sagacity HL SPE product, recoveries for all three metabolites were all above 86%. We tested several different pretreatment method, found out 2:1 serum/IPA ratio could sufficiently release protein bounded vitamin D metabolites, while not enough precipitation to block the SPE. Additionally, certain volume of water needed to obtain 1,25 dihydroxy vitamin D on the Sagacity HL sorbent.

#### Conclusions

Two different types of vitamin D metabolite extraction methods have been established. Vitamin D extraction plate provides reasonable recovery while it is fast and simple. It is very easy to operate by clinical technicians. It is great for 25-hydroxy vitamin D test.

Sagacity HL SPE plate method is high recovery, and still high through put. It could be used for testing low concentration 1,25-dihydroxy vitamin metabolites.