

# Agility

Polymeric Ion Exchange  
Solid Phase Extraction Product



**Agility DVB** ion exchange SPE products are available in four different chemistries: SCX, SAX, WCX, and WAX. Ideal for extraction of a wide range of acidic or basic compounds. Most suited for pharmaceutical, environmental applications and peptide analyses.



## Features

- ▶ Available in varied cartridges and 96-well plate formats
- ▶ Stable over pH 1-14
- ▶ Larger surface area
- ▶ Narrow particle size distribution

## Advantages

- ▶ High capacity for most analytes, allows for use of smaller amount of sorbent, and low elution volume
- ▶ High recoveries and reproducibility
- ▶ Reduced matrix effect



## Applications

- ▶ Pharmaceutical compounds in biological fluids
- ▶ Toxins in foods

