



HPLC Columns & Flash Columns

Orochem Technologies provides a wide range of traditional and customized liquid chromatography columns for reverse phase, normal (polar) phase or ion exchange chemistry. The particle size ranges from 3-200 μm , pore size ranges from 80-300 \AA , and the columns are available in varied sizes



Matrix C18, Monitor C18	Pharmaceuticals	Acetaminophen, Caffeine, 2-Acetamidophenol, Acetanilide, Acetylsalicylic Acid, Azithromycin
Orosil C18	Pharmaceuticals	Tamoxifen, Pregabalin, Budesonide, Eszopiclone, Eprosartan, Stavudine, Ciprofloxacin, Atracurium, Acetaminophen, Caffeine, 2-Acetamidophenol, Acetanilide, Acetylsalicylic Acid
	Vitamins	Riboflavin, Rescorcinol
Prestige C18	Pharmaceuticals	Ampicilline, Chloramphenicol
	Toxins	mycotoxin
	Cosmetics	Sun-screen reagents
	Test Mixtures	Ethyl Benzoate, Toluene, Ethylbenzene, uracil, aniline, p-toluidine, 4-ethyl aniline, N,N-diethylaniline, pyridine, phenol
Reliachrom	Pharmaceuticals, Nucleic Acids	Uracil , 5-fluorouracil, 2-Deoxyuridine, 5-fluoro-2'-deoxyuridine
Orosil Polar C18	Amino Acids	L-Valine, L-Phenylalanine
Reliasil Phenyl		9-Tetrahydrocannabinol
Reliasil C18, C8	Pharmaceuticals, Me-tabolites	Caffeine, Theobromine, Theophylline, Aspirin, Carbamazepine, Aceta-minophen, Amitriptyline, Procainamide
Reliasil C18, C8	Nucleic Acids	Uracil, 5 Flurouracil
	Food Additives/Safety	Benzoic Acid, Melamine
	Organic Acids	Fumaric Acid, Formic Acid, Acetic Acid, Butyric Acid
	Steroids	Estriol, Estradiol, Cortisone, Testosterone, Progesterone, Hydrocortisone
	Vitamins	L-Ascorbic Acid, Nicotnic Acid, Pyridoxine, Nicotinamine, Vitamin D2, Vitamin D3, 25-OH Vitamin D2, 25-OH Vitamin D3, 1,25-Dihydroxy VD3
	Drugs of Abuse, Pharma-ceuticals	Phenobarbital, Secobarbital, Propranolol, MDA, Amphetamine, Methamphetamine, Diazepam, Hydrocodone, Cocaine, Naltrexone, Naloxone, carvedilol, carbamazepine
	Pesticides	Acephate, Dimethoate, Malathion
	Protein and Peptides	Peptide
Reliasil Amine	Alcohols	Isoamyl Alcohol, IPA, Glycerol, Butanol, Ethanol, tagatose, arabinose, galactose, sucrose, rhamnose
Reliasil C30, Matrix-Polychrome	Biofuel, Pharmaceuticals, steroids, fatty acids	Methanol, Ethanol, 2,3-Butanediol, 2-Propanol, carboplatin, decitabine, testosterone & epitestosterone, beta & dexamethasone, oleic acid, palmitic acid, linoleic acid
Reliasil Phenyl-Hexyl	Vitamins	Vitamins E